# **norditalia**<sup>®</sup>

e l e t t r o m e d i c a l i

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## INSTRUCTIONS AND WARNINGS FOR CORRECT AND SAFE USE

**IMPORTANT INSTRUCTIONS.  
KEEP THEM FOR FUTURE  
REFERENCE.**

**THEY SUPPLY IMPORTANT INSTRUCTIONS REGARDING  
OPERATIONS AND SAFETY DURING PRODUCT  
INSTALLATION, USE AND MAINTENANCE.**

### GENERAL WARNINGS

1. After unpacking, make sure that the appliance is intact with no visible signs of damage caused during transport. If in doubt, do not use the appliance and take it back to your retailer.
2. Before connecting the appliance, check that the rating plate specifications, as indicated on the appliance correspond with your mains supply.
3. Do not use adapters, multiple sockets and/or extension cords. If the plug of the appliance is incompatible with the mains socket, you should consult a qualified electrician.
4. This appliance must be used exclusively as a heating pad as indicated in the instructions. Any other use is to be considered hazardous. The manufacturer cannot be held responsible for any damage due to improper or incorrect use.
5. Certain fundamental guidelines must be respected when using any sort of electrical appliance. Specifically:
  - Never touch the appliance with hands or feet that are wet or damp.
  - Never pull on the power cord, or the appliance itself, to remove the plug from the socket.
  - Do not leave the appliance exposed to the elements (rain, sun, etc.).
6. Before cleaning or servicing the appliance, disconnect it by removing the plug from the socket.
7. Should the heating pad or protection fuse break down and/or fail to work properly, turn off the appliance and do not tamper with it. For any repairs, contact your retailer.
8. To avoid any dangerous overheating, we recommend unwinding the power cord completely during use, and disconnecting the plug from the socket when the appliance is not in use.
9. The control unit power cord cannot be replaced. If damaged, it will be necessary to replace the complete unit (chargeable). Contact your retailer.

10. Never immerse the appliance in water.
11. Prolonged use at a high temperature may cause burns. Never fall asleep with the heating pad in operation. Never pull the appliance by the power cord and never use the cord as a grip.
12. This appliance is not designed to be used in hospitals.
13. Do not use the appliance with the infirm, infants or persons insensitive to heat and wearers of pacemakers or other battery-powered medical devices. This appliance is not to be used by young children unless the control has been preset by a parent or guardian or unless the child has been adequately instructed on how to operate the control safely. Children should be supervised to ensure that they do not play with the appliance.
14. Avoid creasing the appliance.
15. This appliance should not be used by children or persons with partially reduced physical, sensorial or mental faculties, and in any case not without specific knowledge and experience, or without the supervision or instructions for use provided by a responsible person.
16. Children must not be allowed to play with the appliance.

**All devices must be disposed of in compliance with the instructions provided by the European Standard n. 2002/96/EC. The recyclable materials used in this appliance can be recovered to prevent environmental degradation. For further information, please contact your local waste disposal authority or your retailer.**

**DO NOT USE IF FOLDED  
OR CRUMPLED**

**DO NOT INSERT  
NEEDLES OR PINS**

# INSTRUCTIONS FOR USE

## IMPORTANT:

This heating pad is fitted with a special safety device which automatically regulates the temperature, by turning the power on and off at regular intervals. The temperature of the heating pad is the maximum allowed by European standards to guarantee user safety.

---

## HOW TO USE

The heating pad is equipped with a control unit with 4 positions:

- 0 off
- 1 minimum temperature
- 2 medium temperature
- 3 maximum temperature

Before turning the heating pad on, check that the control unit is on the 0 position, then put the plug in the socket.

For rapid heating, select level 3; after a few minutes the heating pad will already be pleasantly warm. Place it on the part of the body you wish to heat and select the desired temperature. Do not use the appliance to heat animals.

For best heating results, keep the heating pad covered so that the heat is not lost when the power is turned off.

---

## AUTO OFF FUNCTION

The heating pad is fitted with an automatic shutdown system. After 90 minutes performance time (75 minutes for the New Calor Sabbia and the New Calor Sabbia Cervicale heating pads) the red LED will come on, indicating that the appliance has turned itself off.

To turn the appliance back on, turn the control unit to the "0" position (off) for at least 15 seconds and then re-select the required temperature.

---

## INSTRUCTIONS AND WARNINGS FOR WASHING

### HAND WASH (40°) THE REMOVABLE COVER

The heating pad comes with a separate cover that can be hand washed in warm water (40°) after being removed from the heating pad. We recommend using a mild detergent.

**IMPORTANT:** the appliance itself, cannot be washed or dry cleaned.  
**THE HEATING PAD MUST NOT BE USED IF DAMP.**

---

## STORAGE AND ROUTINE CHECKS

When the appliance is not in use, store it in a dry place and avoid placing heavy objects on it. Avoid creasing the appliance.

Check the appliance and the power cord frequently for any signs of wear or deterioration. Should there be any such signs or if the heating pad fails to work properly, contact your retailer for any possible repairs.

---

## **WARRANTY CONDITIONS**

The warranty is valid for one year from the sale date (two years for European Union countries).

Should the appliance breakdown it should be returned with it's receipt proving the sale date.

The warranty covers the appliance for any manufacturing defects, as long as it hasn't been tampered with.

The warranty does not include:

- breakdown or defects caused by any use other than that foreseen or not complying with the instructions;
- repairs and checks for assumed defects;
- accidental damage such as knocks, falls, etc.

**NOTE:** For repairs and revisions return the appliance in a clean condition.

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